

Pam
+
Nicola

BLURB - re HIV (DIAGNOSTICS) / USE OF PCR

INSTRUCTIONS

- OPEN NCBI
- TYPE ACC # INTO SEARCH
- PULL OUT SECTION OF SEQUENCE WHICH CONTAINS RELEVANT BIT FOR PCR (~500 bases?)

NOW YOU HAVE TO DESIGN PRIMERS FOR THE PCR REACTION TO WORK SUCCESSFULLY.

aatgctgt.....
Sequence

NOTE: NO NUMBERS JUST A STRING OF LETTERS.

OF COURSE DNA IS DOUBLE-STRANDED.

aatgctgt.....
ctacgca.....
double-stranded Sequence

ONCE THIS IS DONE SOMEHOW HIGHLIGHT AREA THAT MUST BE CONTAINED IN THE PCR PRODUCT

- CHOOSE LEFT PRIMER (= MATCH BY CLICK AND DRAG?)

LIST OF 6:
A -
B -
C -
D -
E -
F -

Decide WHICH IS BEST.
(see basic rules of primer design)

YOUR LEFT PRIMER IS INCORRECT! CLUE.....

- CHOOSE RIGHT PRIMER

LIST OF 6:
A —
B —
C —
D —
E —
F —

PRIMERS ARE CHOSEN BUT ARE THEY UNIQUE

- OPEN NCBI
- SEARCH FOR PRIMER SEQUENCE
- HOPEFULLY UNIQUE TO HIV!

END RESULT OF THIS BIT

ds DNA — strands
separate
primers ~~attach~~ bind
new strands created

Basic rules of primer design:

— length — 20-30 bases long
Sequence — avoid long runs of a single letter
eg Actggggggggcgggca

— have c or g at last letter

— avoid self annealing

— must be unique to DNA sequence

GC content — 40-60%

T_m — between 50 and 65°C

$$T_m = 2(A+T) + 4(C+G)$$

