

# User Guide: Primer Design

Ross Barnie  
Dmitrijs Jonins  
Daniel McElroy  
Murray Ross  
Ross Taylor

## Contents

<b>1</b>	<b>About</b>	<b>2</b>
<b>2</b>	<b>Getting Started</b>	<b>2</b>
2.1	System Requirements . . . . .	2
2.2	Downloading and Starting the Application . . . . .	2
2.3	Obtaining a DNA Sequence . . . . .	2
<b>3</b>	<b>The Application</b>	<b>3</b>
3.1	Overview Screen . . . . .	3
3.2	Sequence Entry . . . . .	3
3.3	Target Area Selection . . . . .	3
3.4	Primer Design . . . . .	4
3.5	Melting Temperature . . . . .	4
3.5.1	8th February 2013 Demonstration Only . . . . .	4
<b>4</b>	<b>Known Bugs</b>	<b>4</b>

## List of Figures

1	NCBI Home Page . . . . .	5
2	Searching for a sequence . . . . .	5
3	Search Results . . . . .	5
4	Sequence Page . . . . .	6
5	The Sequence . . . . .	6
6	Overview Screen . . . . .	7
7	DNA Sequence Entry Example . . . . .	7
8	Target Area Selection Example . . . . .	8
9	Primer Selection Example . . . . .	8
10	Primer Design Feedback Example . . . . .	9
11	Melting Temperature Screen Example . . . . .	9

# 1 About

This document will outline how to use the, provisionally named, PrimerDesign application created for use by students of Molecular Methods to learn about Polymerase Chain Reactions (PCR) and primer design.

If you encounter any problems while using the application, feel free to ask your tutor, or email us at [TeamProjectQ@gmail.com](mailto:TeamProjectQ@gmail.com).

## 2 Getting Started

### 2.1 System Requirements

You will need:

- An operating system capable of running Java (you should not need to worry about this as Java runs on almost all modern Windows, Linux and Macintosh operating systems)
- a web browser (such as Google Chrome or Mozilla Firefox)
- Java (at least version 6) installed on your computer (most modern computers have Java, but if not it is freely available online from Oracle)
- Access to the Molecular Methods moodle site by the School of Life Sciences at University of Glasgow
- An internet connection (unless you have downloaded the application and have a DNA sequence ready)

With that we can get started!

### 2.2 Downloading and Starting the Application

First, download the file `PrimerDesign.jar` from the Molecular Methods moodle site and save it somewhere on your computer. Next, double-click on the `PrimerDesign.jar` file you just saved, alternatively if using a Unix-based operating system such as MacOS or Linux you can run the application using the following terminal command:

```
$ java -jar PrimerDesign.jar
```

If you encounter problems here, ensure that you have Java installed on your computer.

### 2.3 Obtaining a DNA Sequence

You will need an internet connection to perform these steps.

Now that you have successfully launched the application, you should prepare a DNA sequence that you wish to manipulate. As an example, we will show you how to obtain the L1CAM sequence, which you should already be familiar with, from the NCBI (National Center for Biotechnology Information) website at <http://www.ncbi.nlm.nih.gov/>.

Use your web browser of choice to go to the NCBI website. You should see something similar to figure 1.

In order to search for a compatible sequence, change the search type to “Nucleotides” from the drop-down menu next to the search bar (highlighted in yellow in figure 2) and search for the sequence you want, for our example this is “L1CAM”, using the search bar (highlighted in red in figure 2).

Now you will be presented with your search results, if you are following our example click on the link highlighted in yellow on figure 3.

Once clicked, you should be presented with something similar to figure 4.

Most of this information is irrelevant to this application, so scroll down until you see the DNA sequence, in our example it should look like figure 5. Now you can simply highlight the sequence (highlighted in yellow in figure 5) and press ctrl-c (or equivalent) to copy the sequence.

### 3 The Application

Please note: the following information and screenshots are subject to change and may not necessarily reflect the current build of the system. Use your best judgement where differences appear.

#### 3.1 Overview Screen

On starting the application (as described in section 2.2), you should see the overview screen (figure 6) with a button on the bottom right labelled “Start”. Once you have finished reading the information on this page, you should press start. Note that you will not be able to return to this page but you will be able to view the Primer Design rules later in the application.

#### 3.2 Sequence Entry

Remember the DNA sequence you copied in section 2.3? Well now is the time to paste it! Simply paste the sequence into the large white area and press the next button. As you can see in figure 7, you do not need to worry about including the “ORIGIN” from the sequence as this will be removed when you press the ‘Next’ button.

#### 3.3 Target Area Selection

Now we have to select what it is we want to copy. To do this, you have to specify the first and the last base of the sequence you wish to copy, using it’s position in the sequence. So if you wish to copy a sequence from position 100 to 500, as in the example on figure 8, you would enter these into the “From” and “To” text boxes.

Note that you can also view the complementary strand, and both strands together, by using the tabs just above where the sequence is.

Note that the sequence is split into seven blocks of 10 bases each (ie 70 bases per line).

### 3.4 Primer Design

You can now see your selected area more clearly and, since primers can include bases from before and after the target, the rest of the sequence is still available to you.

You should design your primers and insert them into the “Forward Primer” and “Reverse Primer” fields, as shown in figure 9, note however that the example data is not designed to be correct.

For the reverse primer, you have a button which will reverse whatever you enter into the “Reverse Primer” field, the “Reverse” button next the field. So if you were to enter `aattccggt`, and press the “Reverse” button, the primer would change to `tggccttaa`.

You can also see the primer design rules again by pressing the “Show Primer Design Rules” button.

When you click next, your primers are checked against the rules described at the start of the application, and you are given a report of where your primers pass and where they fail, if at all. This will look something like 10.

Pressing the “Ok” button will close this window and allow you to continue only if you have passed each rule.

### 3.5 Melting Temperature

Now you are finished! This final screen, which should be similar to figure 11, lets you review your design, showing the melting temperatures of both primers and the primers themselves.

Note that for our example we should not have been allowed to get here due to the feedback we received in figure 10.

#### 3.5.1 8th February 2013 Demonstration Only

Also please note that although we have included the “See Animation” button, it is currently not working so we have disabled it.

## 4 Known Bugs

There are a few known issues with the current build of the system. If you should find any bugs not listed below, please let us know either in person or by e-mail (TeamProjectQ@gmail.com).

- On both the Target Area Selection and Primer Design screens (referred to in sections 3.3 and 3.4) the “Both Strands” tab’s numbers on the left of the screen do not match the lines of the sequence.
- On the Primer Design screen’s (section 3.4) “Both Strands” tab, the highlighting sometimes does not highlight the correct area.




Figure 1: NCBI Home Page




Figure 2: Searching for a sequence




Figure 3: Search Results

**Homo sapiens cell adhesion molecule with homology to L1CAM (close homolog of L1) (CHL1), transcript variant 4, non-coding RNA**

NCBI Reference Sequence: NR\_045572.1

[FASTA](#) [Graphics](#)

[Go to:](#) [Top](#)

LOCUS NR\_045572 2585 bp RNA linear PRI 03-FEB-2013  
DEFINITION Homo sapiens cell adhesion molecule with homology to L1CAM (close homolog of L1) (CHL1), transcript variant 4, non-coding RNA.  
ACCESSION NR\_045572  
VERSION NR\_045572.1 GI:358438320  
KEYWORDS .  
SOURCE Homo sapiens (human)  
ORGANISM [Homo sapiens](#)  
Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Primates; Haplorrhini; Catarrhini; Hominidae; Homo.  
REFERENCE 1 (bases 1 to 2585)  
AUTHORS Long,M.J., Wu,F.X., Li,P., Liu,M., Li,X. and Tang,H.  
TITLE MicroRNA-10a targets CHL1 and promotes cell growth, migration and invasion in human cervical cancer cells  
JOURNAL Cancer Lett. 324 (2), 188-196 (2012)  
PUBMED 22634495  
REMARK GeneRIF: miR-10a expression is upregulated in cervical cancer tissues, and miR-10a promotes cell growth, migration and invasion by targeting CHL1 in human cervical cancer cells.  
REFERENCE 2 (bases 1 to 2585)  
AUTHORS Htt,B., Riordan,S.M., Kukreja,L., Eimer,W.A., Rajapaksha,T.W. and Vassar,R.  
TITLE beta-Site amyloid precursor protein (APP)-cleaving enzyme 1 (BACE1)-deficient mice exhibit a close homolog of L1 (CHL1) loss-of-function phenotype involving axon guidance defects  
JOURNAL J. Biol. Chem. 287 (46), 38408-38425 (2012)  
PUBMED 22989240  
REMARK GeneRIF: BACE1(-/-) axon guidance defects are likely the result of abrogated BACE1 processing of CHL1 and BACE1 deficiency produces a CHL1 loss-of-function phenotype  
REFERENCE 3 (bases 1 to 2585)  
AUTHORS Manning,A.K., Hivert,M.F., Scott,R.A., Grimsby,J.L., Bouatia-Naji,N., Chen,H., Rybin,D., Liu,C.T., Bielak,L.F., Prokopenko,I., Amin,N., Barnes,D., Cadby,G., Hottenga,J.J., Ingelsson,E., Jackson,A.U., Johnson,T., Kanoni,S., Ladenvall,C., Lagou,V., Lahti,J., Lecoeur,C., Liu,Y., Martinez-Larrad,M.T., Montasser,M.E., Navarro,P., Perry,J.R., Rasmussen-Torvik,L.J.,

Figure 4: Sequence Page

[polyA site](#) 2564  
/gene="CHL1"  
/gene\_synonym="CAL1 + L1CAM2"

ORIGIN

```
1 gaaaatgcgc cagcactgga ggggggaatc catggaccgt ggggtgtggt ggttggggag
61 cagcgaatcc cagagagggg ctgagctccc accctcccag accgccgagg ggcgagagg
121 cgcgcggggg cgcgtggttt ccattgctcc ttgtgacgga gccccggggg ctgagatgc
181 cggtcgcgga tggctccctt ttccctcttg tgcctctctc ttctctcgcg ttttttttt
241 ttttttgagt gtgagtggtg gtggggaggg agcacggaga aagtgtatga tttgggaagc
301 agggattgga gccgggaggg tggggaagc cagccctccc gtcccaactc tccggcctcg
361 ctgcgagcca cagtgtctca gccaggggca gggccgggga actcctgcga aaaaccacgg
421 gccggagga cagcagcgcc gggggccag gtgctgaata ctgcaacca taactctgtc
481 ttaatactgc aaacaaatca tagtgaactc aagggaactc taattactg aattctctgt
541 ttggaatttt cctcccgaaa taccaaactc ctcaagaatt taaaacgat ggagtaaaat
601 aatcttcaag ctctgattag attacaacat atgaatgcta attttgttt tacaggtttc
661 aaaccacaaa ttatctgtaa cttagggcac aataagcagt taaactgtta gaaatcttat
721 gagattgagt ttgtgcctat cagatgaata catatttgta tataacattt caattattgc
781 tacatttatg tataatttat ccaaaagcga tgttatccat gcactttaac agaagcctta
841 agtagttgaa tgacagggag tacatttttg ctgcagaaga tcatgataaa acatgtgctg
901 ttgtgtcacc ttggttaagt tggcagactc gttgtatgtt acttatctct acacactgca
961 tctcatactg atttatggat ttggtgactc tattttatgt caaagacata aagagagtag
1021 agaagagaga atagtatttt catttgaatc ttaaaaaagt tctatcctat gagatgttaa
1081 cagctgcctt ctcttgggtg tctgtgaagt tttataaag atattgatta tttccaatg
1141 acacacaaat ctactgtag agaccatcca aatgcttagt tacaccaggg cttgctcctt
1201 gtaccctgag gaggacacga ttaacaaat cctctttttt aaatgttctc aggacaatga
1261 gctgcacata cggcacaaat gcaatttgct tgaccaatct tgcctctgtc tctcttatc
1321 tccgtttaca tgcttcacgc tcagttgaaa ttaattctct acacatttgt ctcatgaaca
1381 catcttttag agcagagttc aagttaaat ttattattca tcaattacct tgagaacaca
1441 atctggggca tattagacgc ttgttaagt tagaggtggg gagaagaag aagaagagaga
1501 gaaacagaag gagggaagaa tttctctgat agatcacttt ctggacactt tctgtctgta
1561 gaagttctga gattgcctct ctggatcctg ctatcttgcc agtcttgcta agagtttcca
1621 tttaactata aacatgttta ctgatgcgtc ttaattgttc cttaactatc aggtttctga
1681 tagcccccac tagccaagtg ggaccctggt caagaagttc acctctgtta aataaatatg
1741 catagatttt caacatctct tcacatgcac tcagaggtga tgaagatttc cactggatgc
1801 tagaacattg tgcataaata atcgaagtgt tatgtctctc agcaatgatg ctaccccata
1861 tctctttttt ttctttccca tgatatgagg atatgagtcg agcagaacta gctccatgcc
1921 ctgtcccgag aattaatcat tctatccatc atttgatgtc ttactctttt aaatatttat
1981 gggctgttat gttcctttct tgtgttccac ccaactctgg gacattcaga agcctcacat
2041 ttcttcttta gccaaacacc ccattccctt catctctctc cagctgtttc ttggtattac
2101 tccaaattgc atacataatt ttggtatgga ggtgccaaaa atgcatttgt ttcaactgga
2161 gttttattga agttgctttt actcaaatca agtacacaaa taaatttgat ttccaaggca
2221 tgaaaattca gttaggctct cgaagggaag ctgaggtata cttatatttc attataatgc
2281 tacaggtaga tttaagtca ctaggaaatc gacactcaca aatcaaaaga ttataatgt
2341 agcactatat gaagaattac attttaagc atagttctgt gatcattttg taagtttttg
2401 ttgttactgt cagtttatgt gtttttatta ataatttat tgaggtataa ctgatacaca
2461 aacgacatat ttaattgtga caatttgatg agtttgaca tatccgagga ctcatgagac
2521 tgtcactgca atcaaggtaa taaatgtatt tattacttcc aaaaaaaaaa aaaaaaaaaa
2581 aaaaaa
```

Figure 5: The Sequence




Figure 6: Overview Screen




Figure 7: DNA Sequence Entry Example

**Primer Design**

### Target Sequence Selection

Specify the region to be amplified by PCR. Indicate the region using the numbers on the left hand side and write these in the boxes below.

From:  To:

	DNA Sequence	Complementary	Double Stranded
1	gaaaatgccg	ccagcctgga	gcggggaatc
71	cagagagggg	ctagagctcc	acgtcccccag
141	ccatggctcc	tgtagcgga	gccggggggg
211	tgctctctc	ttctctcgc	ttttttttt
281	aagtgtataa	tttggggagc	agggattgga
351	tccgggctcg	ctgcgagcca	cagtgtctca
421	gccggaggag	cagcagccgc	gcgggcccag
491	aaacaaatca	tagtgaact	aaggggaact
561	tacccaacct	ctcaagaatt	ttaaaacgat
631	atgaatgcta	attttgtttt	tacaggttcc
701	ttactgttta	aaaatcttat	gagattgagt
771	caattattgc	tacatttatg	tataatttat
841	agtagttgaa	tgacaggagc	tacatttttg
911	ttgggtaagt	tgccagagct	gcttgatggt
981	ttggatgact	tattttatgt	caaaagacata
1051	ttaaaaaagt	tctatcttat	gagatgttaa
1121	atattgatta	ttttccaatg	acacacaaat
1191	cttggtctct	gtaccctctg	gaggacacga
1261	gctgcacata	cggcacaaat	gcaatttgct
1331	aagttaaat	ttattttcca	tcttaacctc
1401	tagaagttgg	gagaaaaaag	aaacaaagaa
1471			

Back Next

Figure 8: Target Area Selection Example

**Primer Design**

### Primer Selection

You now have to choose the forward and reverse primers to amplify the region. Manually type or copy and paste the desired primer sequence into the boxes below. Click the "Show Primer Design Rules" button below to see general primer design rules.

Forward Primer:  Reverse Primer:

	DNA Sequence	Complementary	Double Stranded
1	cttttacggc	ggtcgacct	cgcccttag
71	gtctctcccc	gatctcgagg	tggaggggtc
141	gggtaccgag	accatgcctc	cgggcccccc
211	acggagagag	aaagagagcg	aaaaaaataa
281	ttcactaatt	aaacccctcg	ttccataact
351	agcccgagc	gagctcggtg	gtcacagagt
421	cggtctctcc	gtcgtcgccg	cgccgggttc
491	tttttttatt	ataacattga	ttgccccctg
561	atgggttggg	gagttcttaa	aattttgcta
631	tacttacgat	taaaacaaaa	atgtccaaag
701	aattgacaat	ctttagaata	ctctaaacta
771	gttaataacg	atgttaataa	atatttaata
841	tcatcaactt	actgtccctg	atgttaaaac
911	aaacctattc	acggctctga	cgaactacca
981	aacctactga	ataaaatata	gtttctgtat
1051	aattttttca	agataggata	ctctacaatt
1121	tataactaat	aaaaaggata	tggtgtgta
1191	gaaccgagga	catgggactc	ctctgtgtgt
1261	cgacgtgtat	gcgtgtttta	cggttaacga
1331	acgaagtacg	agtcacattt	aattaagaga
1401	ttcaatttaa	aaataaaagt	agtaaatgag
1471			

Back Show Primer Design Rules Next

Figure 9: Primer Selection Example






Figure 10: Primer Design Feedback Example




Figure 11: Melting Temperature Screen Example